

Exo Iris[®]

Version 3.0 Mobile



INTRODUCTION

Exo Iris[®] is a category-redefining 3-in-1 handheld ultrasound device offering a complete medical imaging ecosystem for use, well, anywhere. Powered by our proprietary silicon technology, Iris integrates advanced imaging, intuitive workflow and real-time artificial intelligence (AI) to deliver immediate answers at the point of care.

WHAT'S NEW

Probe Indicator Orientation

The probe indicator for the cardiac preset is now configurable, defaulting to the right.

WHAT'S FIXED


- The home pinned section now retains user customizations.

KNOWN ISSUES

- For iPad on OS version 26, when the multi-gesture is set to Stage Manager, swiping left to open Mimosa may not work. The workaround is to enable the Mimosa from the control center.

Exo Works[®]

Version 3.0 Web and Mobile



INTRODUCTION

Exo Works[®] is an intuitive POCUS workflow software designed around clinicians. Exo Works[®] helps physicians document and bill exams quickly, and easily manage QA from anywhere.

WHAT'S NEW

- AI-powered web search to quickly find answers to your exam-related questions.
- Export all exams in CSV format (Web only).
- New Exo Works Enterprise license tier:
 - Streamlines billing, educational, and QA workflows.
 - Multi-facilities and specialties support.
 - Create custom exam types and reports.

WHAT'S FIXED

Mobile App:

- When all scans within an exam are deleted by the user, the exam is no longer shown in the Drafts folder. Instead, the user is shown a prompt to confirm the deletion of the entire exam. If confirmed, the exam is deleted.
- When the user opens and deletes the worksheet in a draft exam, the user is no longer stuck. Instead, the user is brought to the previous screen with exam details.

WHAT'S FIXED

Web App:

- When a Specialty Admin renames a Specialty, the previous name no longer shows in the role assignment list.
- If an Admin attempts to assign a role to the user account, the save button is now only enabled when a specialty is added, clearly indicating to the user that a specialty is required when assigning roles.

SmartLink:

- The Smartlink settings page now limits the Worklist AE title field to 16 characters to prevent connection errors.

KNOWN ISSUES

Mobile App:

- When configuring multiple PACS connections in the Essential license tier, multiple PACS may be set as the default. As a workaround, uncheck the PACS that should not be set as the default one.
- When completing an educational exam in the *Order-Based Workflow*, an error is displayed requiring an order assignment. As a workaround, add an order to complete and sign off on the exam.
- For the *Encounter-Based workflow*, when switching from the educational to the clinical category, the user can sign and submit an exam without attestation. As a workaround, remove and add a new worksheet if the exam category needs to be updated.

KNOWN ISSUES

Web App:

- Any user can add or edit worksheets for completed exams (i.e., this is not restricted when it should be).
- For an Order-Based Workflow, if an order is added and then removed from an educational exam after it has been assigned, it is not available in the worklist. As a workaround, the user will have to resend the order.
- When a user downloads a CSV file, the file name is displayed in UTC instead of the user's current local time zone.
- In Pulse, when a Customer Admin is reassigned, they do not receive a welcome email. As a workaround, add the new Customer Admin as a Specialty Admin in Exo Works before updating the Customer Admin in Pulse.

SmartLink:

- Inactive specialties are synced to Smartlink with no indication that they are inactive. As a workaround, the active/inactive status is displayed on the Exo Works Admin page.

To learn more about Exo Iris® and Exo Works®, please visit exo.inc or contact us at exosupport@exo.inc.

Exo Imaging, Inc. ("Exo") has taken care to ensure the accuracy of this document. However, Exo® assumes no liability for errors or omissions and reserves the right to make changes without further notice to any products herein to improve reliability, function, or design. Performance varies depending on system configuration. No product or component can be absolutely secure.

This document contains proprietary information, which is protected by copyright. No license (express or implied, by estoppel or otherwise) to any intellectual property rights is granted by this document. Unauthorized copying of this document, in addition to infringing copyright, might reduce the ability of Exo to provide accurate and current information to users. No part of this document may be reproduced, modified, distributed or otherwise used without prior written permission.

Exo®, Exo Iris®, Exo Works®, and the Exo dots logo are trademarks of Exo Imaging, Inc. in the U.S. and/or other countries.

Copyright 2025, Exo Imaging, Inc.

DOC-0280 Rev. A

Exo Imaging, Inc. 4201 Burton Drive, Santa Clara, CA 95054 USA

Telephone: 1-833-MEET-EXO (1-833-633-8396)