

Exo Iris®

Version 3.2 Mobile



INTRODUCTION

Exo Iris® is a category-redefining 3-in-1 handheld ultrasound device offering a complete medical imaging ecosystem for use anywhere. Powered by our proprietary silicon technology, Iris integrates advanced imaging, intuitive workflow and real-time artificial intelligence (AI) to deliver immediate answers at the point of care.

WHAT'S NEW

Imaging Improvements

- Improved near-field resolution for abdominal imaging.
- Better definition in motion mode.
- Enhanced details in Cardiac imaging for B-mode, Color, and Pulse Wave Doppler modes.

iOS Support

- The Exo Iris app now requires iOS 26. Devices running iOS 18 or earlier are no longer supported.

WHAT'S FIXED


- The midline marker now remains visible on the live image when the depth or the field-of-view is changed.

KNOWN ISSUES

- On iPad, rotating the device while saving a clip, automatically stops the capture. The last few frames might be corrupted. Exo recommends rotating the device after the clip is saved.
- On iPad, the Login button on the Welcome page becomes unresponsive when rotating the device from portrait to landscape orientation. Exo recommends logging in to the app before rotating the device.
- On iPad, when the stage manager feature is enabled, resizing the window may cause some imaging elements to display incorrectly. Exo recommends keeping the window at maximum size while using the app.

Exo Works®

Version 3.2 Web and Mobile



INTRODUCTION

Exo Works® is an intuitive POCUS workflow software designed around clinicians. Exo Works® helps physicians document and bill exams quickly, and easily manage QA from anywhere.

WHAT'S NEW

Exam Unassignment

- Easily unassign yourself from an exam to keep exam ownership accurate.

Personalized View

- View exam details your way. Arrange your columns and save your preferred view

WHAT'S FIXED

- For Order-Based Workflows, exams without an assigned order now display the missing order error directly in both the card and list views, eliminating the need to open the exam to identify the issue.

KNOWN ISSUES

Mobile App:

- When configuring PACS connections in the Essential license tier, multiple PACS may be set as defaults. As a workaround, clear the default selection for any PACS connections that should not be used as the default.
- For Connect and Enterprise license tiers, once a worklist has been configured for a specialty, it no longer appears on the Settings page. However, the worklist connection remains functional, and orders continue to be available for assignment within that specialty.

KNOWN ISSUES - CONTINUED

Web App:

- When assigning specialty-based roles using the "Bulk Edit" action in User Management, the "All" option for Exam Types is not selected by default. If access to all exam types is intended, manually select "All" before saving.
- An error may occur when removing the Specialty Admin role from a user. As a workaround, first assign the user the Facility Admin role, then remove the Specialty Admin role.

SmartLink:

- SmartLink does not currently display whether a specialty is active or inactive. Administrators can view each specialty's status on the Exo Works Admin page.

To learn more about Exo Iris® and Exo Works®, please visit exo.inc or contact us at exosupport@exo.inc.

Exo Imaging, Inc. ("Exo") has taken care to ensure the accuracy of this document. However, Exo® assumes no liability for errors or omissions and reserves the right to make changes without further notice to any products herein to improve reliability, function, or design. Performance varies depending on system configuration. No product or component can be absolutely secure.

This document contains proprietary information, which is protected by copyright. No license (express or implied, by estoppel or otherwise) to any intellectual property rights is granted by this document. Unauthorized copying of this document, in addition to infringing copyright, might reduce the ability of Exo to provide accurate and current information to users. No part of this document may be reproduced, modified, distributed or otherwise used without prior written permission.

Exo®, Exo Iris®, Exo Works®, and the Exo dots logo are trademarks of Exo Imaging, Inc. in the U.S. and/or other countries.

Copyright 2026, Exo Imaging, Inc.

DOC-0323 Rev. A

Exo Imaging, Inc. 4201 Burton Drive, Santa Clara, CA 95054 USA

Telephone: 1-833-MEET-EXO (1-833-633-8396)